

## Upper Limb

Double Shoulder Brace		
0581	0582	0583
0585	0587	0588

### Description

An orthosis designed to fit around both shoulders, providing compression and support. The two-part construction allows for individual patient anatomy. The brace maintains the head of the humerus in the glenoid fossa. The Neoprene material provides warmth to the shoulder joints and for patient comfort.

### Intended Purpose

To provide upper limb support for an array of injuries and conditions, including Pre or Postoperative treatment for: Lax capsule, Shoulder instability, Rhomboid and Trapezius strains

### Sizing

The size is determined by the upper arm circumference. Measure the circumference at the largest part of the upper arm.

### ⚠ Fitting Instructions

We recommend the initial fitting of the device be conducted by a suitably qualified Healthcare Professional who will advise fitting, removal and period of use.

1. Place one arm into the cuff making sure the hook and loop is on the back of the arm and secure the hook and loop fastening **See Fig.1.**
2. Repeat step one with the opposite arm **See Fig.2.**
3. Secure across the back of the neck and secure the hook and loop fastening **See Fig.3.**
4. To achieve more abduction undo the hook and loop on the back, retighten until the desired effect and re-secure the hook and loop **See Fig.4.**



### **Contraindications**

- Must not be used by individuals for whom compression is contraindicated
- Do not use on open wounds
- Any known allergy and/or hypersensitivity to any material listed in the product composition
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy

### **Warnings and Precautions**

- Carefully read all instructions and warnings before use
- Follow all instructions to ensure proper performance of the brace
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area
- Do not reuse for another patient, doing so risks cross-infection and can compromise product integrity
- To ensure the straps are not over tightened, perform regular skin and circulation checks, especially for patients with diabetes, vascular deficiencies and neurological conditions
- Should any adverse reactions occur, please remove, do not use and contact your healthcare professional or provider
- Check brace and all components for wear and tear
- The hook and loop fasteners must always be fastened when the device is not being worn or when being washed
- The durability of the device may be compromised by certain factors, e.g. objects with sharp edges or damage to the hook and loop fasteners

### **Washing and care instructions**

- Cover all the hook and loop straps
- Hand wash the brace in cool water using a mild detergent
- Rinse and air dry flat
- Do not spin or tumble dry
- Do not use fabric conditioner at any time
- Check for any degradation that may affect use
- Periodically check the hook and loop straps are operating securely

### **Product composition**

Neoprene, Hook and Loop

### **Storage and transport conditions**

Store in a cool, dry place out of direct sunlight, in the original packaging.

### **Recycling and disposal**

Packaging and constituent parts should be recycled or safely disposed of in accordance with local or national laws.

### **Serious incident**

Report any serious incident to the manufacturer and the competent authority of the EU Member State, or Country in which you reside.

### **Wash symbols**



Hand wash



Dry flat



Do not tumble dry



Do not bleach



Do not dry clean



Do not iron



**Beagle Orthopaedic**

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

e: info@beagleortho.com t: +44 (0) 1254 268 788 w: beagleortho.co.uk



This Medical Device conforms to the requirements of Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017



**International Associates Auditing & Certification Limited**

The Black Church, St. Mary's Place, Dublin, D07 P4AX