

# Bea Pad

## INDICATIONS FOR USE

- » Tennis elbow

## MEASUREMENT

- » Forearm circumference

## FITTING INSTRUCTIONS

- » Open the hook & loop fastening and slip arm in to just below the elbow.
- » The round pressure pad, inside the brace, should lie below the bony part on the outside of the elbow, leaving enough room to allow the elbow to bend.
- » Put the strap through the D-ring pull back and attach to fastening.

## WASHING INSTRUCTIONS

- » Hand wash
- » Air dry

### Sizing Information

<b>Size</b>	<b>Forearm Circ (mm)</b>	<b>Code</b>
Small	225 - 275	BTE/03
Medium	275 - 310	BTE/05
Large	310 - 340	BTE/07
XLarge	340 - 380	BTE/08



## Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

E: [info@beagleortho.com](mailto:info@beagleortho.com) W: [beagleortho.com](http://beagleortho.com) T: 01254 268 788 F: 01254 268789



European Healthcare & Device Solutions Ltd  
Sutton House, Bishopstown Road  
Cork, T12 Y9TC, Ireland



Any product supplied by Beagle Orthopaedic is intended for single patient use only. Re-use of any product will invalidate the CE mark and any warranty. In case of any adverse reaction, please cease use immediately and contact your clinician for further advice.