

Bea Wrap/Bea Wrap Long

INDICATIONS FOR USE

- » Increasing medial/lateral support
- » Rheumatoid arthritis
- » Mild to moderate hyperextension
- » Mild to moderate valgus/varus deformity

MEASUREMENT

- » Circumference of knee with leg straight

FITTING INSTRUCTIONS

- » Open all the hook and loop straps
- » Place the brace behind the knee making sure the hole lies at the back of the knee and the hinges lie along the inside and outside of the knee. If required the hinges can be removed and shaped to the leg to provide a snug fit
- » Wrap the brace around the leg and close with the hook and loop pads making sure the knee sits in the space on the front
- » Pass the straps through the D-rings and attach back on themselves to tighten, the straps may need to be tightened when standing to prevent slippage

*Hinges should not be bent within 25mm of centre of hinge

WASHING INSTRUCTIONS

- » Remove metal hinges
- » Hand wash
- » Air dry



Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR
E: info@beagleortho.com W: beagleortho.com T: 01254 268 788 F: 01254 268789



European Healthcare & Device Solutions Ltd
Sutton House, Bishopstown Road
Cork, T12 Y9TC, Ireland



Any product supplied by Beagle Orthopaedic is intended for single patient use only. Re-use of any product will invalidate the CE mark and any warranty. In case of any adverse reaction, please cease use immediately and contact your clinician for further advice.