

Paediatric Wrist Brace

INDICATIONS FOR USE

- » Healing Fractures
- » Strains and sprains
- » Juvenile RA

MEASUREMENT

- » Wrist circumference at narrowest point

FITTING INSTRUCTIONS

- » Undo all hook & loop straps
- » Remove aluminium palmar bar and bend to required shape. Place bar along underside of wrist onto palm. If there are spaces bend as required
- » Reinsert palmar bar and cover with end tab
- » Place brace on underside of wrist and align palmar stay to palm of hand
- » Fasten wrist hook & loop straps to ensure snug fit

WASHING INSTRUCTIONS

- » Remove metal bar
- » Hand wash
- » Air dry



Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR
E: info@beagleortho.com W: beagleortho.com T: 01254 268 788 F: 01254 268789



European Healthcare & Device Solutions Ltd
Sutton House, Bishopstown Road
Cork, T12 Y9TC, Ireland



Any product supplied by Beagle Orthopaedic is intended for single patient use only. Re-use of any product will invalidate the CE mark and any warranty. In case of any adverse reaction, please cease use immediately and contact your clinician for further advice.