

Knee & Lower Limb

Bea Lok Brace			
0713	0715	0717	0718

Description

An orthosis featuring push-lock knee hinges to provide automatic locking at full extension, simple push-button allows for ease in everyday movements e.g sitting. The brace also features a front opening design for ease of donning with 4 circumferential straps. The standard length is 420mm, with a shell lining in plush material that provides padding for extra comfort.

Intended Purpose

To provide knee and lower limb support for an array of injuries and conditions, including Pre or Postoperative treatment for: Trauma, Rigid knee support for ambulation, CVA or other neurological impairments, Knee stability for degenerative conditions

Sizing

Measure the circumference of the knee in full extension.

⚠ Fitting Instructions

We recommend the initial fitting of the device be conducted by a suitably qualified Healthcare Professional who will advise fitting, removal and period of use.

1. Unfasten all the straps and open the brace and lay the brace flat **See Fig.1.**
2. Offer the brace up and ensure alignment is maintained **See Fig.2.**
3. Secure the four straps in any order (tip - in order to limit migration we suggest that the straps closest to the knee are secured first), further adjustments/tightening can be made when the patient is standing if appropriate **See Fig. 3.**
4. The knee hinges are simple push locks, use your finger to lock and unlock, locking provides automatic full extension, unlocking to allow for ease of everyday movement such as sitting **See Fig. 4.**



Contraindications

- Must not be used by individuals for whom compression is contraindicated
- Do not use on open wounds
- Any known allergy and/or hypersensitivity to any material listed in the product composition
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy

Warnings and Precautions

- Carefully read all instructions and warnings before use
- Follow all instructions to ensure proper performance of the brace
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area
- Do not re-use for another patient, doing so risks cross-infection and can compromise product integrity
- To ensure the straps are not over tightened, perform regular skin and circulation checks, especially for patients with diabetes, vascular deficiencies and neurological conditions
- Should any adverse reactions occur, please remove, do not use and contact your healthcare professional or provider
- Check brace and all components for wear and tear
- The hook and loop fasteners must always be fastened when the device is not being worn or when being washed
- The durability of the device may be compromised by certain factors, e.g. objects with sharp edges or damage to the hook and loop fasteners
- Do not bend the hinges within 25mm of the hinge centre

Washing and care instructions

- Cover all the hook and loop straps and remove metal hinges
- Hand wash the brace in cool water using a mild detergent
- Rinse and air dry flat
- Do not spin or tumble dry
- Do not use fabric conditioner at any time
- Check for any degradation that may affect use
- Periodically check the hook and loop straps are operating securely
- Before using the brace again, replace the metal hinges

Product composition

Bioprene, Aluminium, Hook and Loop, Plastic Buckle

Storage and transport conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and disposal

Packaging and constituent parts should be recycled or safely disposed of in accordance with local or national laws.

Serious incident

Report any serious incident to the manufacturer and the competent authority of the EU Member State, or Country in which you reside.

Wash symbols



Hand wash



Dry flat



Do not tumble dry



Do not bleach



Do not dry clean



Do not iron




Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

e: info@beagleortho.com t: +44 (0) 1254 268 788 w: beagleortho.co.uk



This Medical Device conforms to the requirements of Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017



International Associates Auditing & Certification Limited

The Black Church, St. Mary's Place, Dublin, D07 P4AX