

## Knee & Lower Limb

Tri Panel Economy Immobiliser			
ETPI12/BB	ETPI14/BB	ETPI16/BB	ETPI18/BB
ETPI20/BB	ETPI22/BB	ETPI24/BB	ETPI12/BG
ETPI14/BG	ETPI16/BG	ETPI18/BG	ETPI20/BG
ETPI22/BG	ETPI24/BG		

### Description

An orthosis made of a durable three piece construction that allows for easy application to a wide variety of leg lengths. The brace contains independent stays on each panel to enable shaping for a more customised fit. The immobiliser features four circumferential straps to aid in prevention of migration and positioning straps to aid in the fitting procedure.

### Intended Purpose

To provide knee and lower limb support for an array of injuries and conditions, including Pre or Postoperative treatment for: Early cast removal, Patella dislocation, Sprains and strains of knee joint, ACL, PCL, MCL, LCL injuries, Fractures, Arthritis

### Sizing

One size fits most, the support is universal in size and varied in length.

### △ Fitting Instructions

We recommend the initial fitting of the device be conducted by a suitably qualified Healthcare Professional who will advise fitting, removal and period of use.

1. Undo all hook and loop straps and remove the side steels. The side steels should also be shaped to accommodate the patient's anatomy. **See Fig.1.**
2. The side panels can be detached and re-positioned to accommodate the patient's leg circumference. **See Fig.2.**
3. Close the panel straps to secure the brace to the patient's leg. **See Fig.3.**
4. Wrap the circumferential straps around the brace and secure. The straps may need to be adjusted once the patient is standing. **See Fig.4.**



### **Contraindications**

- Must not be used by individuals for whom compression is contraindicated
- Do not use on open wounds
- Any known allergy and/or hypersensitivity to any material listed in the product composition
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy

### **Warnings and Precautions**

- Carefully read all instructions and warnings before use
- Follow all instructions to ensure proper performance of the brace
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area
- Do not reuse for another patient, doing so risks cross-infection and can compromise product integrity
- To ensure the straps are not over tightened, perform regular skin and circulation checks, especially for patients with diabetes, vascular deficiencies and neurological conditions
- Should any adverse reactions occur, please remove, do not use and contact your healthcare professional or provider
- Check brace and all components for wear and tear
- The hook and loop fasteners must always be fastened when the device is not being worn or when being washed
- The durability of the device may be compromised by certain factors, e.g. objects with sharp edges or damage to the hook and loop fasteners

### **Washing and care instructions**

- Cover all the hook and loop straps and remove metal bars
- Hand wash the brace in cool water using a mild detergent
- Rinse and air dry flat
- Do not spin or tumble dry
- Do not use fabric conditioner at any time
- Check for any degradation that may affect use
- Periodically check the hook and loop straps are operating securely
- Before using the brace again, replace the metal bars

### **Product composition**

8mm Velfoam with Towelling Lining, Aluminium, Hook and Loop

### **Storage and transport conditions**

Store in a cool, dry place out of direct sunlight, in the original packaging.

### **Recycling and disposal**

Packaging and constituent parts should be recycled or safely disposed of in accordance with local or national laws.

### **Serious incident**

Report any serious incident to the manufacturer and the competent authority of the EU Member State, or Country in which you reside.

### **Wash symbols**



Hand wash



Dry flat



Do not tumble dry



Do not bleach



Do not dry clean



Do not iron




**Beagle Orthopaedic**

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

e: info@beagleortho.com t: +44 (0) 1254 268 788 w: beagleortho.co.uk



This Medical Device conforms to the requirements of Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017



**International Associates Auditing & Certification Limited**

The Black Church, St. Mary's Place, Dublin, D07 P4AX