

Spinal

Wiltshire Maternity Pelvic Belt

WMB/096

WMB/116

WMB/130

WMB/146

Description

An orthosis to support throughout pregnancy. The Wiltshire Maternity Pelvic Belt was designed as a slimline maternity belt for comfort when sitting. The belt provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly. The brace features an elastic back section with two elastic side straps for variable tension.

Intended Purpose

To provide back and pelvic support for an array of injuries and conditions, including Pre or Postoperative treatment for: Pubic Support, Sacroiliac support throughout all stages of pregnancy.

Sizing

Measure the circumference of the hips.

⚠ Fitting Instructions

We recommend the initial fitting of the device be conducted by a suitably qualified Healthcare Professional who will advise fitting, removal and period of use.

1. Undo the hook and loop fastening **See Fig.1.**
2. Align the belt to the midline of the lower back, with the flexible steels lying on either side of the spine **See Fig.2.**
3. Wrap the soft elastic underneath the belly alleviating the weight on the lower spine and lifting the belly and close the hook and loop fastening **See Fig.3.**
4. The side straps provide additional support, pull the straps out and across **See Fig.4.**
5. Secure underneath the belly using the hook and loop fastening **See Fig.5.**

Fig.1.



Fig.2.



Fig.3.



Fig.4.



Fig.5.



Contraindications

- Must not be used by individuals for whom compression is contraindicated
- Do not use on open wounds
- Any known allergy and/or hypersensitivity to any material listed in the product composition
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy

Warnings and Precautions

- Carefully read all instructions and warnings before use
- Follow all instructions to ensure proper performance of the brace
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area
- Do not reuse for another patient, doing so risks cross-infection and can compromise product integrity
- To ensure the straps are not over tightened, perform regular skin and circulation checks, especially for patients with diabetes, vascular deficiencies and neurological conditions
- Should any adverse reactions occur, please remove, do not use and contact your healthcare professional or provider
- Check brace and all components for wear and tear
- The hook and loop fasteners must always be fastened when the device is not being worn or when being washed
- The durability of the device may be compromised by certain factors, e.g. objects with sharp edges or damage to the hook and loop fasteners

Washing and care instructions

- Cover all the hook and loop straps and remove metal bars
- Hand wash the brace in cool water using a mild detergent
- Rinse and air dry flat
- Do not spin or tumble dry
- Do not use fabric conditioner at any time
- Check for any degradation that may affect use
- Periodically check the hook and loop straps are operating securely

Product composition

Velfoam, Elastic, Cotton, Hook and Loop, Latex

Storage and transport conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.


Recycling and disposal

Packaging and constituent parts should be recycled or safely disposed of in accordance with local or national laws.

Serious incident

Report any serious incident to the manufacturer and the competent authority of the EU Member State, or Country in which you reside.


Wash symbols




Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

e: info@beagleortho.com **t:** +44 (0) 1254 268 788 **w:** beagleortho.co.uk

 This Medical Device conforms to the requirements of Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017

 **International Associates Auditing & Certification Limited**
The Black Church, St. Mary's Place, Dublin, D07 P4AX