

Knee & Lower Limb

Paediatric Tri Panel		
TP08.GRN	TP08.RED	TP08.BLU
TP10.GRN	TP10.RED	TP10.BLU

Description

An orthosis made from an 8mm Velfoam with grey towelling lining and available in 8" and 10" lengths. The knee brace provides effective support and stabilisation to an affected area, helping to reduce pain and aid recovery. The durable three piece construction allows for easy application to a wide variety of leg sizes. The independent stays on each panel enable a more customised fit.

Intended Purpose

To provide knee and lower limb support for an array of injuries and conditions, including Pre or Postoperative treatment for: Early cast removal, Patella dislocation, Sprains and strains of knee joint, ACL, PCL, MCL, LCL injuries, Fractures, Arthritis

Sizing

One size fits most, the support is universal in size and varied in length.

⚠ Fitting Instructions

We recommend the initial fitting of the device be conducted by a suitably qualified Healthcare Professional who will advise fitting, removal and period of use.

1. Undo all hook and loop straps and remove the side steels. The side steels should also be shaped to accommodate the patient's anatomy. **See Fig.1.**
2. The side panels can be detached and re-positioned to accommodate the patient's leg circumference. **See Fig.2.**
3. Close the panel straps to secure the brace to the patient's leg. **See Fig.3.**
4. Wrap the circumferential straps around the brace and secure. The straps may need to be adjusted once the patient is standing. **See Fig.4.**



Contraindications

- Must not be used by individuals for whom compression is contraindicated
- Do not use on open wounds
- Any known allergy and/or hypersensitivity to any material listed in the product composition
- Caution should be used if prescribed to patients with diabetes, vascular deficiency and neuropathy

Warnings and Precautions

- Carefully read all instructions and warnings before use
- Follow all instructions to ensure proper performance of the brace
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area
- Do not reuse for another patient, doing so risks cross-infection and can compromise product integrity
- To ensure the straps are not over tightened, perform regular skin and circulation checks, especially for patients with diabetes, vascular deficiencies and neurological conditions
- Should any adverse reactions occur, please remove, do not use and contact your healthcare professional or provider
- Check brace and all components for wear and tear
- The hook and loop fasteners must always be fastened when the device is not being worn or when being washed
- The durability of the device may be compromised by certain factors, e.g. objects with sharp edges or damage to the hook and loop fasteners

Washing and care instructions

- Cover all the hook and loop straps and remove metal bars
- Hand wash the brace in cool water using a mild detergent
- Rinse and air dry flat
- Do not spin or tumble dry
- Do not use fabric conditioner at any time
- Check for any degradation that may affect use
- Periodically check the hook and loop straps are operating securely
- Before using the brace again, replace the metal bars

Product composition

Velfoam with Towelling Lining, Aluminium, Hook and Loop

Storage and transport conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and disposal

Packaging and constituent parts should be recycled or safely disposed of in accordance with local or national laws.

Serious incident

Report any serious incident to the manufacturer and the competent authority of the EU Member State, or Country in which you reside.

Wash symbols



Hand wash



Dry flat



Do not tumble dry



Do not bleach



Do not dry clean



Do not iron




Beagle Orthopaedic

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR

e: info@beagleortho.com t: +44 (0) 1254 268 788 w: beagleortho.co.uk



This Medical Device conforms to the requirements of Regulation (EU) 2017/745 of the European Parliament and of the council of 5 April 2017



International Associates Auditing & Certification Limited

The Black Church, St. Mary's Place, Dublin, D07 P4AX