

CTLSO Airback Combo - Specification Sheet

A customised option to provide maximum support for high thoracic control

standard female sizes

sizing information					
size	chest circ (mm)	waist circ (mm)	hip circ	tick to select	model No
X-small	660 - 760	550 - 660	780 - 890	<input type="checkbox"/>	B/2099
small	680 - 780	600 - 710	830 - 930	<input type="checkbox"/>	B/2100
medium	730 - 830	660 - 760	890 - 990	<input type="checkbox"/>	B/2101
large	840 - 940	760 - 860	1010 - 1115	<input type="checkbox"/>	B/2102
X-large	910 - 1010	915 - 1015	1060 - 1170	<input type="checkbox"/>	B/2103

standard male sizes

sizing information					
size	chest circ (mm)	waist circ (mm)	hip circ	tick to select	model no
X-small	840 - 940	710 - 810	840 - 940	<input type="checkbox"/>	B/2104
small	860 - 940	710 - 810	830 - 940	<input type="checkbox"/>	B/2105
medium	890 - 995	840 - 940	965 - 1065	<input type="checkbox"/>	B/2106
large	990 - 1090	910 - 1060	1015 - 1115	<input type="checkbox"/>	B/2107
X-large	1060 - 1165	1040 - 1140	1115 - 1220	<input type="checkbox"/>	B/2108



sizing information		
size	chin to shoulder depth (mm)	tick to select
short	57	<input type="checkbox"/>
regular	76	<input type="checkbox"/>
tall	96	<input type="checkbox"/>
extra tall	114	<input type="checkbox"/>

Please provide chin to waist measurement

.....