

Bea Flex Wrist Brace

Material - Foam Laminate



Manufactured
in the UK



Designed in
Association with
Clinicians



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. Made from a plush foam laminate material and designed in collaboration with Occupational Therapists. A rigid support featuring the Dual Y shaped strapping which creates a pulley effect to provide excellent compression around the whole wrist for maximum support and comfort. With removable metal palmer bar.



Indications for use

- » Tendinitis and tendon disorders of the wrist
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis

Features and benefits

- » Rigid support featuring the Dual Y shaped strapping which creates a pulley effect to provide excellent compression
- » Maximum support and comfort
- » Can be hand washed and air dried

Sizing Information

Size	Wrist Circ (mm)	Code Left	Code Right
Small	130 - 150	BEAF/33L	BEAF/33R
Medium	150 - 170	BEAF/35L	BEAF/35R
Large	180 - 200	BEAF/37L	BEAF/37R
X-Large	200 - 230	BEAF/38L	BEAF/38R

Beagle
ORTHOPAEDIC