

BEA Wrap

Material - Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. The Bea Wrap is designed with front opening contrasting proximal and distal strap to counter rotation. It has removable 10° hyperextension stop hinges. The open popliteal area provides added comfort during flexion and the easy grip tabs aid patients with poor dexterity. Universal: fits left or right.



Indications for use

- » Increasing medial/lateral support
- » Rheumatoid arthritis
- » Mild to moderate hyperextension
- » Mild to moderate valgus/varus deformity

Features and benefits

- » Front opening
- » Universal: left or right
- » Contrasting proximal and distal strap to counter rotation
- » Removable 10° hyperextension stop
- » Open popliteal area for added comfort during flexion
- » Easy grip tabs for patients with poor dexterity
- » Can be hand washed and air dried

Sizing Information

Size	Knee Circ	Code
Small	270 - 355	BWST/03
Medium	330 - 405	BWST/05
Large	380 - 455	BWST/07
XLarge	430 - 500	BWST/08
XXLarge	500 - 555	BWST/09
XXXLarge	550 - 605	BWST/10
XXXXLarge	600 - 655	BWST/11