

Darco Medsurg Healing Shoe



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Darco MedSurg Shoe has a semi-rigid sole that is built around a MetaShank forefoot support that provides control under the metatarsal heads and protection after osteotomies. The rocker sole reduces plantar pressure on the forefoot and heel by over 25%. The square toe design acts as a bumper to provide additional protection when K-wires are present. A strapless fastening system can be removed and rotated to provide a left or right overlap and the new ankle strap eliminates heel slippage. A zoned outsole provides traction under the midfoot. A removable insole is included but for a more customised fit the peg assist insole is available.



Indications for use

- » Used postoperatively for stability and off-loading, e.g. after foot osteotomy procedures

Features and benefits

- » Breathable mesh upper
- » Removable insole
- » Can be worn on the left or right foot

Sizing Information - Men

Size	Shoe Size	Code
Small	up to 8	MQM1B
Medium	8.5 - 10	MQM2B
Large	10.5 - 12	MQM3B
XLarge	12.5 - 14	MQM4B

Sizing Information - Women

Size	Shoe Size	Code
Small	Up to 4	MQW1B
Medium	4.5 - 5.5	MQW2B
Large	6 - 7.5	MQW3B