

Darco WCS[®] Light Wound Care Shoe



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

The Wound Care Shoe System Light from DARCO offers the proven four-layer-insole system with individual customisation for offloading included in a closed shoe. The ultra soft & breathable upper provides maximum comfort and leaves extra space for bandages. Featuring a seamless toe area with extra space for light bandages, with an elastic & breathable upper and comfortable fitting provided by high heel cup. Sold in pairs and 5 sizes available in black (adults).



Indications for use

- » Treatment for diabetic foot syndrome
- » Healing shoe for open wounds and mild ulcerations
- » Designed for pressure redistribution
- » Used for postoperative phase after toe and/or forefoot surgery

Features and benefits

- » Includes customisable insoles
- » Replacement “four-layer” insole available
- » Sold in pairs

Sizing Information

Size	Shoe Size	Code
XSmall	Up to 3	WCSL0
Small	4 - 5	WCSL1
Medium	7.5 - 8.5	WCSL2
Large	9 - 10	WCSL3
XLarge	11 - 12	WCSL4