

Flexiform Wrist Brace

Material - Ventilated Elastic



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

Manufactured in the UK. Made in a breathable elastic material that provides support to the whole wrist region for added reassurance. Maintains full finger movement for rehabilitation exercises. Easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons.



Indications for use

- » Healing Fractures
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis
- » Carpal Tunnel syndrome

Features and benefits

- » Webbed elastic material
- » Provided support to the whole wrist region
- » Relieves strain on tendons
- » Helps to improve circulation and promote healing
- » Reduces inflammation to aid with tingling and numbness
- » Can be hand washed and air dried

Sizing Information

Colour - Beige

Size	Wrist Circ (mm)	Code Left	Code Right
Small	150 - 170	FFEB/3L.BE	FFEB/3R.BE
Medium	170 - 180	FFEB/5L.BE	FFEB/5R.BE
Large	180 - 200	FFEB/7L.BE	FFEB/7R.BE
XLarge	200 - 230	FFEB/8L.BE	FFEB/8R.BE



Sizing Information

Colour - Black

Size	Wrist Circ (mm)	Code Left	Code Right
Small	150 - 170	FFEB/3L.BLK	FFEB/3R.BLK
Medium	170 - 180	FFEB/5L.BLK	FFEB/5R.BLK
Large	180 - 200	FFEB/7L.BLK	FFEB/7R.BLK
XLarge	200 - 230	FFEB/8L.BLK	FFEB/8R.BLK