

Foam Stirrup Ankle

Material - Polypropylene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

The Foam Stirrup Ankle is comprised of 2 polypropylene outer shells connected with a soft, adjustable foot piece. The unique inner memory foam pads provide uniform compression and a customised fit. The pads are wider at the bottom to prevent biting. It provides stability to the ankle by minimising ankle rotation. Universally sized, one size fits all.



Indications for use

- » Chronic ankle instabilities and prophylactic use

Features and benefits

- » Provides stability to the ankle by minimizing ankle rotation
- » Comprises 2 polypropylene outer shells connected with a soft, adjustable foot piece
- » Unique memory foam pads provide uniform compression and a customised fit
- » Pads are wider at the bottom to prevent 'biting'
- » Universal, one size fits all

Sizing Information

Size	Code
Universal	US8813F