

# Hallux Valgus Splint

Material - Bioprene



Manufactured  
in the UK



Anti Bacterial



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

Manufactured in the UK, from BIOPRENE a new generation of antibacterial material treated with AEGIS TM Microbe Shield; which gives long-lasting protection against bacterial build up, making it ideal when the product is worn continuously for extended periods. The breathable BIOPRENE material provides compression for full contact support. The hook and loop straps make it very easy to don and doff and the abduction strap for the first metatarsal joint (Great Toe) supports the toe alignment. Contains a removable aluminium support stay for corrective positioning. Can be worn Post-Operatively or as a Night Resting Splint.



## Indications for use

Provides superior post operative toe alignment after following procedures:

- » Hallux Valgus
- » Hammer Toe
- » Tailor's Bunion

## Features and benefits

- » Bioprene is a new generation of antibacterial material treated with AEGIS TM microbe shield
- » Aegis gives long-lasting protection against bacterial build up, ideal when the product may be worn continuously for extended periods
- » Breathable Bioprene provides compression for full contact support
- » Easy fasten Hook & Loop straps
- » Abduction strap for the first metatarsal joint (great toe) for toe alignment
- » Aluminium support stays for corrective positioning
- » Worn post-operatively or as a Night Resting Splint

### Sizing Information

Size	Male Shoe Size	Female Shoe Size	Left	Right
Small	-	2 - 4	HVS/3L	HVS/3R
Medium	-	4 - 6	HVS/5L	HVS/5R
Large	5 - 7	6 - 7	HVS/7L	HVS/7R
XLarge	8 - 11	8 - 11	HVS/8L	HVS/8R