

Hinged Knee Sleeve

Material - Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. Open patella design to keep patella in position. It provides compression and support for weak and injured knees and retains body heat to stimulate blood flow and accelerate healing. The reinforced patella stabiliser to prevent displacement, with removable 10° hyperextension stop hinges. Universal: will fit left or right.



Indications for use

- » Patella disorders that need efficient control

Features and benefits

- » Open patella support relieves pressure on the knee cap
- » Provides compression and support for weak and injured knees
- » Retains body heat to stimulate blood flow and accelerate healing
- » Reinforce patella stabiliser to prevent displacement
- » Removeable 10° hyperextension stop hinges
- » Universal: left or right
- » Can be hand washed and air dried

Sizing Information

Size	Thigh Circ	Knee Circ	Calf Circ	Code
Small	380 - 420	310 - 340	300 - 330	HK/223
Medium	400 - 440	330 - 370	320 - 360	HK/225
Large	440 - 500	370 - 420	360 - 400	HK/227
XLarge	460 - 520	400 - 460	400 - 440	HK/228