

Neoprene Thumb Spica

Material - Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. Made in neoprene that allows for gentle compression to the whole wrist region effectively reducing pain and inflammation. The metal stay provides support and protection to the thumb which is maintained in a neutral position.



Indications for use

- » CMC ligament sprains and strains
- » RSI
- » Arthritis

Features and benefits

- » Gentle compression to the whole wrist region for added reassurance
- » Protection to the thumb which is maintained in a neutral position
- » Targeted firm support to relieve pain
- » Can be hand washed and air dried

Sizing Information

Size	Wrist Circ (mm)	Code Left	Code Right
Small	130 - 170	NTS/03L	NTS/03R
Medium	170 - 210	NTS/05L	NTS/05R
Large	210 - 250	NTS/07L	NTS/07R
X-Large	250 - 290	NTS/08L	NTS/08R

Beagle
ORTHOPAEDIC