

Standard Ankle Sleeve

Material - Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK, from a Neoprene material, the standard ankle sleeve provides compression, comfort and support to the ankle injury. The hook and loop fastening system makes it very easy to don and doff.



Indications for use

- » Acute ankle injuries
- » Arthritis

Features and benefits

- » Compression, gives both comfort & support to the ankle injury
- » Incorporates hook & loop closure for easier application

Sizing Information

Size	Circ of ankle (mm)	Code
Small	180 - 200	SAS/03
Medium	200 - 220	SAS/05
Large	220 - 240	SAS/07
XLarge	240 - 270	SAS/08