

Bea Shoulder Sling

Material - Soft Cotton



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. The Bea shoulder sling provides immobilisation following injury or surgery with an integrated neck pad provides additional patient comfort. Featuring an adjustable neck and waist strap. Elbow mobilisation is possible without complete sling removal and the thumb loop effectively maintains sling placement. The soft cotton material provides increased patient comfort.



Indications for use

- » Immobilisation following injury or surgery

Features and benefits

- » Integrated adjustable neck pad provides additional patient comfort
- » Adjustable neck and wrist strap
- » Elbow mobilisation is possible without complete sling removal
- » Thumb loop effectively maintains sling placement
- » Soft cotton material for increased patient comfort
- » Can be hand washed and air dried



Sizing Information

Size	Length of Sling (mm)	Codes
Small	340 - 380	0663
Medium	370 - 430	0665
Large	420 - 480	0667
X-Large	480 +	0668