

# Orthomerica California Hip System



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Comfortable, total-contact reinforced mesh pelvic component based on the popular California Compression Orthosis supports the back and aligns pelvis and hip joint. Universal pelvic, thigh and hip joint components reduce inventory and make stocking all sizes more cost efficient. Hip joint settings provide a strong flexion/ extension stop with mild, adjustable resistance to adduction. Comfortable mesh thigh component. Patient-adjustable compression to increase comfort during activities of daily living. Abduction wedge regulates the clearance of the hip joint over the trochanter and incision site. Lightweight, low profile components are easier to wear under clothing which can increase patient compliance and wearing time.



## Indications for use

- » Post-op revision patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals
- » Adult patients with CVA, spinal cord injuries or head trauma can benefit from the California Hip with bilateral thigh cuffs to widen the base of support in standing and sitting
- » Patients who have been non-compliant with rigid plastic hip orthoses
- » Patients whose cognitive function impairs their ability to follow hip precautions
- » Patients with weak hip abductors who need a kinesthetic cue and positioning to prevent dislocation
- » Elderly patients confined to bed who are at risk for dislocation

# Product Codes

## California Hip Thigh Component

### Sizing Information

Size	Prox Thigh (mm)	Distal Circ (mm)	Code
Small	380 - 530	305 - 432	B/3976
Medium	430 - 610	330 - 450	B/3977
Large	480 - 680	350 - 500	B/3978
XLarge	530 - 730	430 - 580	B/3979

## California Hip Pelvic Component

### Sizing Information

Size	Waist Circ (mm)	Hip Circ (mm)	Code
Small	560 - 710	660 - 810	B/3970
Medium	710 - 860	810 - 960	B/3971
Large	860 - 1010	960 - 1110	B/3972
XLarge	1010 - 1160	1110 - 1270	B/3973

## California Hip Hinge

### Sizing Information

Size	Code
V5 extended	B/4327
Short	B/4326*

\*V5 Short hinge - Patients between 4'8" & 5'2"

## Please State Left or Right Configuration

### Sizing Information

Config	Code
Left	CHA/L
Right	CHA/R

## Features and benefits

- » Comfortable, total-contact reinforced mesh pelvic component based on the popular California™ Compression Orthosis supports the back and aligns pelvis and hip joint
- » Universal pelvic, thigh and hip joint components reduce inventory and make stocking all size more cost efficient
- » Hip joint settings provide a strong flexion/extension stop with mild, adjustable resistance to adduction
- » Comfortable mesh thigh component
- » Patient-adjustable compression to increase comfort during activities of daily living
- » Abduction wedge regulates the clearance of the hip joint over the trochanter and incision site
- » Lightweight, low profile components are easier to wear under clothing which can increase patient compliance and wearing time
- » Can be hand washed and air dried