

Clasby Shoulder Humeral Brace

Material - Bioprene



Manufactured in the UK




Designed in Association with Clinicians



Hand Wash & Air Dry Only



Same Day Dispatch



Breathable



ISO:13485 Accredited

Manufactured in the UK. Clinically proven, the clasby humeral brace is effective at controlling pain and is associated with high patient satisfaction scores. Constructed from a Bioprene which is a new generation of anti bacterial material treated with AEGIS microbe shield. It provides compression for full contact support. Left and right in design. Available with metal or plastic stays.



Indications for use

- » Mid and proximal diaphyseal fractures

Features and benefits

- » Maintains integrity of the skin more than traditional rigid braces
- » The chest strap distributes weight across the body whilst not impinging on the breast area, especially in females
- » It supplies traction and provides mild abduction to control and manage pain
- » Easy to apply, enhancing compliance and promote early rehabilitation
- » There are 6 easy to remove metal stays that exert pressure equally to stabilise the fracture
- » Trimmable around the elbow to provide a more customised fit
- » There are 2 dual pull compression straps to apply compression and prevent rotation of the fracture
- » Extension piece provided to accommodate oedema
- » Can be hand washed and air dried

Sizing Information

Metal Stays		Codes	
Size	AC Joint to Elbow (mm)	Left	Right
X-small	280	CBS/2L	CBS/2R
Small	300	CBS/3L	CBS/3R
Medium	320	CBS/5L	CBS/5R
Large	340	CBS/7L	CBS/7R
X-large	350	CBS/8L	CBS/8R

Sizing Information

Plastic Stays* none stock item, please allow 5 days for delivery		Codes	
Size		Left	Right
XSmall		CBS/32L	CBS/32R
Small		CBS/33L	CBS/33R
Medium		CBS/35L	CBS/35R
Large		CBS/37L	CBS/37R
XLarge		CBS/38L	CBS/38R

Beagle
ORTHOPAEDIC