

Clavicle Splint

Material - Back Panel: White Laminate

Straps: Webbing



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. This splint helps maintain ideal posture by pulling evenly over the scapulas. Unnoticeable when worn under clothing and the slip resistant buckles ensure patient compliance. To order, measure from centre of shoulder blades, over shoulder, around axilla & back to centre.




Indications for use

- » Clavicle fracture
- » Clavicle injury
- » General posture support
- » Shoulder retraction

Features and benefits

- » Helps maintain ideal posture by pulling evenly on shoulder
- » Unnoticeable when worn under clothing
- » Slip resistant buckle ensure patient compliance
- » Measure from centre of shoulder blades, over shoulder around axilla & back to centre



Sizing Information

Size	Shoulder - Axilla (mm)	Codes
Small	560 - 640	CLA/03
Medium	740 - 720	CLA/05
Large	720 - 800	CLA/07
X-Large	800 - 880	CLA/08