

Darco PegAssist[®] Insole

For Medsurg, Relief Dual, OrthoWedge & HeelWedge



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

The perfect choice for the treatment of diabetic ulcerations because it minimises pressure on wounds for faster healing while allowing patients to maintain their mobility. Also eliminates the lever effect on osteotomies of the MT1 by offloading the first MPJ. The insole can also be cut flush with the metatarsal heads to float the toes when K-Wires are present. Provides short-term pain management for patients suffering due to bony abnormalities (e.g. prominent metatarsal heads) or other painful conditions that may later be surgically corrected. The pegs can be removed to off-load pressure from any area of the plantar surface of the foot. A single-use marking tube is included to mark the area of desired pressure reduction and to ensure proper peg removal. The peg support board adheres to black side of insole to keep pegs from collapsing and the off-loaded area is kept intact. For use with the Medsurg Healing Shoe, Relief Dual, OrthoWedge & HeelWedge.



Indications for use

Selective off-loading for:

- » Plantar lesions
- » Plantar ulcerations
- » Pressure points

Features and benefits

- » Insole is designed for use with the right or left foot

Sizing Information - Men Medsurg, Relief Dual

Size	Shoe Size	Code
Small	Up to 8	PTQM1
Medium	8.5 - 10	PTQM2
Large	10.5 - 12	PTQM3
XLarge	12.5 - 14	PTQM4

Sizing Information - Universal for Men & Women OrthoWedge & HeelWedge

Size	Shoe Size	Code
XSmall	Up to 8	PQ0
Small	8.5 - 10	PQ1
Medium	10.5 - 12	PQ2
Large	12.5 - 14	PQ3
XLarge		PQ4

Sizing Information - Women Medsurg, Relief Dual

Size	Shoe Size	Code
Small	Up to 4	PTQW1
Medium	4.5 - 5.5	PTQW2
Large	6 - 7.5	PTQW3

Beagle
ORTHOPAEDIC