

Darco Relief Dual® Off-Loading Shoe



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

The Relief Dual® innovations revolutionise the off-loading technology. This newly designed Dual sole combines the benefits of a EVA shock absorbing material and the anti-slip features of rubber; for a ground-breaking step in the post-operative treatment for foot care. In addition, through these improvements Darco have optimised the rolling characteristics of the shoe. Breaking from convention, they have increased the surface area of the rigid shank to the entire sole area; maximising pressure re-distribution while effectively immobilising the metatarsophalangeal joint I & V (MTP) and the distal interphalangeal joint I & V (DIP).



Indications for use

- » Used postoperatively for stability and off-loading e.g. after foot osteotomy procedures

Features and benefits

- » Innovative dual-sole technology
- » Optimal impact relief - slip-resistant rubber sole - light weight breathable upper
- » Rigid shank throughout entire innersole
- » Wide opening for bulky dressings
- » Easy on/off, and easy insole access
- » Can be worn on the left or right foot
- » Low profile - no height compensation necessary

Sizing Information - Men

Size	Shoe Size	Code
Small	Up to 8	RD-M1
Medium	8.5 - 10	RD-M2
Large	10.5 - 12	RD-M3
XLarge	12.5 - 14	RD-M4

Sizing Information - Women

Size	Shoe Size	Code
Small	Up to 4	RD-W1
Medium	4.5 - 5.5	RD-W2
Large	6 - 7.5	RD-W3