

Double Shoulder Brace

Material - Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Manufactured in the UK. Made from neoprene to provide warmth to shoulder joints and for patient comfort. The two part construction allows for individual patient anatomy. This brace maintains the head of the humerus in the glenoid fossa.



Indications for use

- » Shoulder instability
- » Lax capsule
- » Rhomboid and Trapezius Strains

Features and benefits

- » Two part construction to allow for individual patient anatomy
- » Maintains head of humerus in the glenoid fossa
- » Offers compression and support to both shoulders
- » Neoprene construction for patient comfort
- » Provides warmth to the shoulder joint
- » Fits under clothing
- » Can be hand washed and air dried



Sizing Information

Size	Upper arm circ (mm)	Codes
Small	280 - 330	0583
Medium	330 - 380	0585
Large	380 - 430	0587
X-Large	430 - 480	0588