

Flexiform Wrist Thumb Brace

Material - Ventilated Elastic/Neoprene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

Manufactured in the UK. A ventilated webbed elastic material. Elastic provides support to the whole wrist region for added reassurance. Wrist brace with palmar bar and thumb extension bar which can be shaped to provide a more custom fit. Provides additional support to the first CMC and MCP joints and the removable stay can be shaped.



Indications for use

- » Healing Fractures
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis

Features and benefits

- » Provided support to the whole wrist region
- » Relieves strain on tendons
- » Ventilated elastic body to allow air to circulate
- » Neoprene thumb piece to provide full compression for comfort
- » Can be hand washed and air dried

Sizing Information

Colour - Beige

Size	Wrist Circ (mm)	Code Left	Code Right
Small	150 - 170	FFWT/3L.BE	FFWT/3R.BE
Medium	170 - 180	FFWT/5L.BE	FFWT/5R.BE
Large	180 - 200	FFWT/7L.BE	FFWT/7R.BE
XLarge	200 - 230	FFWT/8L.BE	FFWT/8R.BE



Sizing Information

Colour - Black

Size	Wrist Circ (mm)	Code Left	Code Right
Small	150 - 170	FFWT/3L.BLK	FFWT/3R.BLK
Medium	170 - 180	FFWT/5L.BLK	FFWT/5R.BLK
Large	180 - 200	FFWT/7L.BLK	FFWT/7R.BLK
XLarge	200 - 230	FFWT/8L.BLK	FFWT/8R.BLK