

# Humeral Extended Fracture Brace

Material - LD Polyethylene



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Extended functional brace indicated for closed diaphyseal fractures. The shoulder cap design limits distal brace migration giving greater protection and compliance. Easily trimmed for optimum fit and to allow for range of motion in the elbow. With adjustable hook & loop straps to maintain desired compression.



## Indications for use

- » Closed proximal and mid shaft diaphyseal humeral fractures

## Features and benefits

- » Extended functional brace indicated for closed diaphyseal fractures
- » Shoulder cap design limits distal brace migration giving greater protection and compliance
- » Easily trimmed for optimum fit and to allow for range of motion in the elbow
- » Adjustable hook & loop straps maintain desired compression
- » Universal left and right
- » Item includes: 1 fracture brace, 2 fracture socks



### Sizing Information

Size	Proximal Circ (mm)	Distal Circ (mm)	Code
XSmall	200 - 250	150 - 200	B/1089
Small	230 - 280	200 - 250	B/1090
Medium	280 - 330	250 - 280	B/1091
Large	330 - 380	280 - 330	B/1092
XLarge	380 - 430	330 - 380	B/1093