

Orthomerica Newport 4



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Low-profile design to improve patient comfort and compliance with a padded posterior section for enhanced comfort and improved lordotic fit. No proximal/medial plastic on thigh component to interfere with hygiene. The thigh component design adjusts with clamping disc to accommodate varying shapes and to prevent knee flexion impingement. Preset, quick release closures ensure same fit each time. A patented clamping technology on pelvic band adjusts to male or female anatomy and a wide range of circumferences. The universal pelvic design fits right and left to reduce inventory. Also available with universal style thigh cuff.



Indications for use

- » Post-op hip revision patients
- » Primary arthroplasty patients at risk to dislocate, e.g. patients with congenital hip dysplasia or spastic/tight adductors
- » As a prophylaxis to reinforce hip precautions
- » Inoperable patients at risk

Features and benefits

- » Patented offset adjustment limits contact with incision site
- » Low profile design to improve patient comfort and compliance
- » Padded posterior section for enhanced comfort and compliance
- » No proximal/medial plastic on thigh component to interfere with hygiene
- » Thigh component design adjusts with clamping disk to accommodate varying shapes and sizes
- » Preset, quick release closures ensure same fit each time
- » Patented clamping technology on pelvic band adjusts to male or female anatomy and a wide range of circumferences
- » Universal pelvic design fits right and left to reduce inventory
- » Thigh extends down to medial part of knee to encourage abduction
- » Can be hand washed and air dried

Product Codes

Newport® 4 Thigh Component

Sizing Information

Size	Prox Thigh (mm)	Distal Circ (mm)	Left	Right
Small	380 - 530	300 - 430	B/3770.06	B/3775.06
Medium	430 - 610	330 - 460	B/3771.06	B/3776.06
Large	480 - 680	350 - 510	B/3772.06	B/3777.06
XLarge	530 - 730	430 - 580	B/3773.06	B/3778.06

Newport® 4 Pelvic Component

Sizing Information

Size	Waist Circ (mm)	Hip Circ (mm)	Code	Replacement Pads
Small	560 - 710	660 - 810	B/3870	B/3812.02
Medium	710 - 860	810 - 960	B/3871	B/3811.02
Large	860 - 1010	960 - 1110	B/3872	B/3810.02
XLarge	1010 - 1160	1110 - 1270	B/3873	B/3813.02

Newport® 4 Hinges

Sizing Information

Size	Code
V5 Short curved hinge	B/4320*
V5 Curved hinged	B/4323

*V5 Short hinge - Patients between 4'8" & 5'2"

Thigh Replacement Liners

Sizing Information

Size	Left	Right
Small	B/3820.02	B/3825.02
Medium	B/3821.02	B/3826.02
Large	B/3822.02	B/3827.02
XLarge	B/3823.02	B/3828.02