

# Russell Hall Sling

Material - Velfoam



Manufactured  
in the UK



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Manufactured in the UK. This high arm sling can maintain the arm at 45 or 90 degrees. The waist strap limits shoulder movement and carries weight of the arm across the back and shoulder for enhanced comfort. Universal: one size will fit on left or right arm.



## Indications for use

- » Post trauma
- » Post operative
- » Oedema control
- » Paralysis
- » CVA

## Features and benefits

- » Easy to don and doff
- » Universal: one size fits all
- » Universal: left and right
- » Maintains arm at 45° or 90°
- » Waist strap limits shoulder movement & carries the weight of the arm across the shoulder for enhanced comfort
- » Colour coded tabs to aid in fitting the product once patient/carer is in own environment
- » Can be hand washed and air dried

### Sizing Information

Size	Code
Universal	0305/RH