

# Sporlastic Swede-O Ankle Lok Brace

Material - Latex Free



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Exclusive ANKLE LOK® offset panel traps the laces between the inner and outer layer to maintain the position for long periods of time. The close spaced eyelets concentrate holding power where it's needed most. The full elastic back panel encourages blood flow to the Achilles' tendon and virtually eliminates the chance for blistering. With its internal U-shaped spiral stays it provides extra support and further minimise the chance for ankle injury. While the Arch fits the contour of the foot and is seamless so it virtually eliminates irritation to the bottom of the foot.



## Indications for use

- » Conservative therapy for sprained ankles
- » Post-operative follow-up treatment

## Features and benefits

- » Exclusive Ankle Lok® offset panel traps the laces between the inner and outer flap to hold the laces together longer than any other brace
- » Exclusive close spaced eyelets concentrate the holding power where its needed most creating the most effective heel lock
- » Full elastic back ensures complete unrestricted blood flow to the Achilles' tendon and virtually eliminates the chance for blistering
- » Internal U-shaped spiral stays provide extra support and further minimise the chance for ankle injury
- » Arch fits the contour of the foot and is seamless so it virtually eliminates irritation to the bottom of the foot

### Sizing Information

Size	Ankle Circ (cm)	Code
2	20 - 22	070722
3	22 - 24	070723
4	24 - 26	070724
5	26 - 28	070725
6	28 - 31	070726