

Orthomerica Fixed Ankle Orthosis



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Standard AFO design pre-fabricated from polypropylene to a 1/8" thickness for necessary rigidity and comfort. Indicated for flaccid or spastic foot drop. Can be heated for further modification and an individual fit. A very versatile design containing posterior ribs for increased support.



Indications for use

- » Flaccid foot drop with spasticity

Features and benefits

- » Prefabricated from polypropylene to a 1/8" thickness for necessary rigidity and comfort
- » Indicated for post-stroke flaccid or spastic foot drop
- » Can be heated for further modification and individual fit
- » Very versatile design
- » Posterior ribs for increased support
- » Velcro® brand strap with roller buckle to optimise fit

Sizing Information

Size	Shoe Size	Footplate Length (mm)	Brace Height (mm)	Left	Right
Small	5 - 7	160	300	B/3540	B/3545
Medium	7.5 - 10	160	340	B/3541	B/3546
Large	10.5 - 13	170	340	B/3542	B/3547
XLarge	13+	190	355	B/3543	B/3548