

Paediatric Hip Abduction Splint Adjustable



Hand Wash &
Air Dry Only



ISO:13485
Accredited

Lightweight paediatric orthosis offers advantage of abduction adjustment over the standard paediatric hip splint model. Sturdy polypropylene frame set in 90 degree hip flexion - abduction easily adjustable 40 degree through 180 degree by aluminium slotted screw assembly. Fits infants from 3 months [6-12 lbs.] to 18 months [over 30 lbs.]. X-large size effective for certain cerebral palsy patients. It has more total contact than pavlik harness. Easy application, removal facilitates nappy changes and the closed-cell foam does not retain moisture.



Indications for use

- » Hip Dysplasia
- » XL effective for some cerebral palsy patients

Features and benefits

Lightweight paediatric orthosis offers advantage of abduction adjustment over the standard paediatric hip splint model. Sturdy polypropylene frame set in 90 degree hip flexion - abduction easily adjustable 40 degree through 180 degree by aluminium slotted screw assembly

- » Fits infants from 3 months (6-12lbs) to 18 months (over 30 lbs). XLarge size effective for certain cerebral palsy patients.
- » More total contact than Pavlik harness
- » Easy application, removal facilitates nappy changes
- » Closed-cell foam does not retain moisture

Sizing Information

Size	Avg. Age	Weight guide (IBS)	Code
XSmall	Up to 3 mths	6 - 12	B/3906
Small	3 - 6 mths	13 - 16	B/3907
Medium	6 - 12 mths	17 - 21	B/3908
Large	1 year	22 - 30	B/3909
XLarge	18 mths	30 +	B/3910

NOTE: Adjustable version is a non-stock item and is non-returnable

Beagle
ORTHOPAEDIC