

# Orthomerica Supra-Lite™ AFO



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop. Injection moulded polypropylene allows for different thicknesses throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility. It can be heat modified for minor fitting adjustments. The dorsum strap is included. It has removable / washable soft laminate calf padding for improved comfort and hygiene.



## Indications for use

- » Flaccid foot drop
- » Post stroke patients

## Features and benefits

- » Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop
- » Injection moulded polypropylene allows for different thickness throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility
- » Can be heat modified for minor fitting adjustments
- » Dorsum strap included
- » Removable/washable soft laminate calf padding for improved comfort and hygiene
- » Can be hand washed and air dried only

### Sizing Information

Size	Footplate Length (mm)	Brace Height (mm)	Code Left	Code Right
Small	240	265	B/3500	B/3505
Medium	280	285	B/3501	B/3506
Large	280	330	B/3502	B/3507
XLarge	300	390	B/3503	B/3508