

# Soft Collar

Material - Foam



Manufactured  
in the UK



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

Manufactured in the UK. Ideal for whiplash pain relief, neck strain and arthritis of the neck available in a range of lengths and depths to obtain the perfect fit. Adjustable fasteners maintain the desired compression to obtain some much-needed pain relief.

To order please **measure the depth from under the ear to the collar bone** and **measure the circumference around the neck and add 4 inches** for length measurement.



## Indications for use

- » Neck strain/sprain
- » Arthritis
- » Whiplash injuries

## Features and benefits

- » Adjustable fasteners maintain the desired compression to obtain some much-needed pain relief
- » Maintain head position providing stabilisation for comfort
- » Helps to control pain by restricting position

### Sizing Information

Size (Inches)	Code
17" X 2.5"	SFC/171
17" X 3"	SFC/172
17" X 3.5"	SFC/173
18" X 2.5"	SFC/181
18" X 3"	SFC/182
18" X 3.5"	SFC/183
19" X 2.5"	SFC/191
19" X 3"	SFC/192
19" X 3.5"	SFC/193
20" X 2.5"	SFC/201
20" X 3"	SFC/202
20" X 3.5"	SFC/203
22" X 2.5"	SFC/221
22" X 3"	SFC/222
22" X 3.5"	SFC/223

# Soft Collar Black

## Sizing Information

Size (Inches)	Code
17" X 2.5"	SFC/171.BLK
17" X 3"	SFC/172.BLK
17" X 3.5"	SFC/173.BLK
18" X 2.5"	SFC/181.BLK
18" X 3"	SFC/182.BLK
18" X 3.5"	SFC/183.BLK
19" X 2.5"	SFC/191.BLK
19" X 3"	SFC/192.BLK
19" X 3.5"	SFC/193.BLK
20" X 2.5"	SFC/201.BLK
20" X 3"	SFC/202.BLK
20" X 3.5"	SFC/203.BLK
22" X 2.5"	SFC/221.BLK
22" X 3"	SFC/222.BLK
22" X 3.5"	SFC/223.BLK

