

Wiltshire Maternity Belt

Material -



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

Manufactured in the UK. The slimline design of this belt makes it comfortable to wear when sitting. It provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly, featuring an elastic back section with 2 elastic side straps, for variable tension.



Indications for use

- » Sacroiliac support throughout all stages of pregnancy

Features and benefits

- » Pelvic support belt with an elastic back section
- » 2 elastic side straps, for variable tension
- » The slimline design of this belt makes it comfortable to wear when sitting
- » Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly
- » Can be hand washed and air dried

Sizing Information

Hip Circ (mm)	Code
960 - 1160	WMB/096
1160 - 1130	WMB/116
1300 - 1460	WMB/130
1460 - 1600	WMB/146