

Paediatric D-Ring Wrist Brace

Material - Cotton




Manufactured in the UK



Hand Wash & Air Dry Only



Same Day Dispatch



Breathable



ISO:13485 Accredited

Manufactured in the UK. The paediatric D-ring wrist brace saves clinical time as no cast removal required. A contoured design to improve fit with 3 touch and close fastenings. The stockinette tongue to ensure optimal comfort. Available in a variety of colours and patterns to suit any requirements.



Indications for use

- » Green stick fracture
- » Buckle fracture
- » Healing wrist fractures
- » Healing ligament injuries
- » Post injury

Features and benefits

- » Saves clinical time as no cast removal required
- » Contoured design to improve fit
- » 3 touch and close fastenings
- » Stockinette tongue to ensure optimal comfort
- » Variety of colours and patterns to suit any requirements
- » Can be hand washed and air dried

Sizing Information - Blue

Size	Length (mm)	Wrist Circ (mm)	Code Right	Code Left
Small	130	130	DR/3RB	DR/3LB
Medium	150	140	DR/5RB	DR/5LB
Large	175	155	DR/7RB	DR/7LB
XLarge	180	185	DR/8RB	DR/8LB

Sizing Information - Red

Size	Length (mm)	Wrist Circ (mm)	Code Right	Code Left
Small	130	130	DR/3RR	DR/3LR
Medium	150	140	DR/5RR	DR/5LR
Large	175	155	DR/7RR	DR/7LR
XLarge	180	185	DR/8RR	DR/8LR

Beagle
ORTHOPAEDIC

Sizing Information - Purple

Size	Length (mm)	Wrist Circ (mm)	Code Right	Code Left
Small	130	130	DR/3RP	DR/3LP
Medium	150	140	DR/5RP	DR/5LP
Large	175	155	DR/7RP	DR/7LP
XLarge	180	185	DR/8RP	DR/8LP

Sizing Information - Green

Size	Length (mm)	Wrist Circ (mm)	Code Right	Code Left
Small	130	130	DR/3RG	DR/3LG
Medium	150	140	DR/5RG	DR/5LG
Large	175	155	DR/7RG	DR/7LG
XLarge	180	185	DR/8RG	DR/8LG



Custom D-ring Brace

For a custom made D-ring wrist brace, simply choose your required size and state colours.

Sizing Information

Size	Length (mm)	Wrist Circ (mm)	Code Right	Code Left
Small	125	120	DR/3R	DR/3L
Medium	150	140	DR/5R	DR/5L
Large	180	160	DR/7R	DR/7L
XLarge	180	180	DR/8R	DR/8L

Colours

Material	Hook & Loop	Binding
Black	Black	Black
White	White	White
Red	Red	Red
Green	Green	Green
Blue	Blue	Blue
Yellow		
Purple		