

Bea Cool Wrist Brace

Material - Bioprene



Manufactured
in the UK




Anti-Bacterial



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

Manufactured in the UK. Made from a ventilated Bioprene that allows moisture to wick away from the skin and contains anti-bacterial properties to eliminate odours. Black material may be preferred for workplace use. The removable easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons.



Indications for use

- » Healing fractures
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis

Features and benefits

- » Ventilated bioprene to allow moisture to wick away from the skin
- » Black material may be preferred for workplace use
- » Anti-bacterial properties to eliminate odours
- » Easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons
- » Can be hand washed and air dried

Sizing Information

Size	Wrist Circ (mm)	Code Left	Code Right
Small	150 - 170	BWB/23L	BWB/23R
Medium	170 - 180	BWB/25L	BWB/25R
Large	180 - 200	BWB/27L	BWB/27R
X-Large	200 - 230	BWB/28L	BWB/28R