

Orthomerica Paediatric UFO



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Introducing an improved prefabricated night splint for the effective treatment of Plantar Fasciitis and Achilles Tendinitis. The Paediatric UFO™ comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Load Prolonged Stretch [LLPS] to the soft tissues of the plantar fascia and Achilles tendon. LLPS promotes gentle stretching, healing and re-modelling of these tissues, effectively reducing the stress and pain caused by inflammation and tissue shortening. Such support can help break the cycle of injury, pain and re-injury due to otherwise unlimited contraction and deformity of the soft tissues at night. The simple variable dorsiflexion design allows the patient to quickly and easily adjust the amount of applied stretch to provide the most effective pain relief. Adjustable bilateral dorsiflexion assist straps with a CoolFoam™ replaceable liner that wicks moisture away from the skin for enhanced hygiene. Universal will fit left or right. Clinical papers are available upon request.



Indications for use

- » Plantar Fasciitis
- » Achilles Tendonitis
- » Sever's disease
- » Idiopathic toe walking
- » Cerebral palsy & other Neuromuscular disease

Features and benefits

- » The UFO™ comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Loas Prolonged Stretch (LLPS) to the soft tissues of the plantar fascia and Achilles tendon
- » LLPS promotes gentle stretching, healing and re-modelling of these tissues, effectively reducing the stress and pain
- » Such support can help break the cycle of injury, pain and re-injury due to otherwise unlimited contraction and deformity of the soft tissue at night
- » The simple variable dorsiflexion design allows the patient to quickly and easily adjust the amount of applied stretch
- » For pattern options see custom Fabrication swatch

Sizing Information

Size	Footplate Length (mm)	Brace Height (mm)	Code	Colour
Small	140	165	B/3614	White
Medium	170	250	B/3615	White
Large	190	265	B/3616	White
Medium	170	250	B/3615.02	Butterfly
Large	190	265	B/3616.02	Butterfly
Medium	170	250	B/3615.03	Tornado
Large	190	265	B/3616.03	Tornado
Custom	-	-	B/310/PE	White
Custom	-	-	B/310/PE/T	Pattern

Beagle
ORTHOPAEDIC