

# Sporlastic Achillodyn


Material - Latex Free



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

The Achillodyn provides proprioceptive muscular stabilisation. Constructed from a special fabric that has high absorption properties without the floating effect as with full silicone. Contains two heel wedges one to temporarily reduce the tension and the other for the healthy side to avoid uneven pelvis. The separate silicone pad gently massages the oedema of the paratenon during movement and increases the blood circulation.



## Indications for use

- » Acute, sub-acute and chronic achillodynia
- » Achilles tendinitis

## Features and benefits

- » Separate heel wedge to temporarily reduce tension
- » Additional heel wedge for the healthy side to avoid scoliotic pelvis
- » Special fabric with high absorption properties without the floating effect as with full silicone
- » Separate silicone pad gently massages the oedema of paratenon during movement and increases blood circulation
- » Proprioceptive muscular stabilisation

### Sizing Information

Size	Circ of leg above ankle (mm)	Code
1	up to 190	0770711
2	190 - 210	0770712
3	210 - 230	0770713
4	230 - 250	0770714
5	250 - 270	0770715