

# Standard Wrist Brace 7"

Material - Tri Weave




Manufactured  
in the UK



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

Manufactured in the UK and made of a tri-weave elastic. The elastic provides support to the whole wrist region for added reassurance. Full finger movement is maintained for rehabilitation exercises. Removable easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons. 4 strap fastening system. Sold singular or boxes of 6.



## Indications for use

- » Healing Fractures
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis

## Features and benefits

- » Elastic Provides support to the whole wrist region for added reassurance
- » Tri-weave contains latex but not in contact with the skin
- » Maintains full finger movement for rehabilitation exercises
- » Easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons
- » 4 finger tabs
- » Provides support to reduce strain on tendons
- » Can be hand washed and air dried

### Sizing Information

Size	Wrist Circ (mm)	Code Box of 6		Code Singles	
		Left	Right	Left	Right
Small	150 - 170	SWB/3L	SWB/3R	SWB/3L/SGL	SWB/3R/SGL
Medium	170 - 180	SWB/5L	SWB/5R	SWB/5L/SGL	SWB/5R/SGL
Large	180 - 200	SWB/7L	SWB/7R	SWB/7L/SGL	SWB/7R/SGL
X-Large	200 - 230	SWB/8L	SWB/8R	SWB/8L/SGL	SWB/8R/SGL

**Beagle**  
ORTHOPAEDIC