

BioSkin Crossfire Conforma Pull On


Material - Sirius



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Crossfire Conforma is a supportive and lightweight hinged knee brace which is ideal for people recovering from a variety of knee injuries who want to remain mobile and active while they recover. An extremely lightweight and flexible hinge offers increased joint stability while allowing for a full range of motion. The base sleeve of the CrossFire Conforma pull on provides compression to help reduce swelling, provide pain relief and enhance proprioception. The CrossFire Conforma pull on is also a great knee brace for sports and injury prevention.



Indications for use

- » Common injuries
- » ACL & PCL injury
- » MCL & LCL injury
- » Meniscus tears
- » Knee instability
- » Hyperextension
- » Knee sprains

Features and benefits

- » Features lightweight and compressive material over the patella as shown in the wrap around versions donning sleeve
- » These braces also feature a more generous cut above the knee to accommodate a wider variety of leg shapes

Sizing Information

Circumference with knee fully extended

Size	Mid-patella (mm)	180mm above mid - patella (mm)	Code
Small	330 - 360	460 - 510	61861
Medium	360 - 410	510 - 560	61862
Large	410 - 430	560 - 640	61863
XLarge	430 - 460	640 - 710	61864

Contact customer service for other sizes.

Beagle
ORTHOPAEDIC