

# BioSkin DP2


**Material** - Soft Microfleece



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

This premium wrist brace is a comfortable way to stabilise the wrist joint while recovering from carpal tunnel syndrome or wrist sprains. It is recommended to wear a wrist brace at night to keep the wrist in a neutral position while sleeping, which is why BioSkin uses soft, moisture-wicking microfleece around the skin and across the sensitive area between your thumb and fingers to keep you cool, dry, and comfortable. The lightweight, supportive stays on the top and bottom of the wrist are removable and can be reshaped for a customizable fit. Constructed with straps instead of laces, means the brace can easily be put on with one hand. The low-profile and functional design of the brace allows for full use of your fingers while limiting the motion of your wrist.



## Indications for use

- » Tendinitis and tendon disorders of the wrist
- » Sprains
- » Repetitive strain injuries
- » Rheumatoid arthritis

## Features and benefits

- » Customisable aluminium stays can be bent for comfort and specific positioning or removed if not needed
- » Soft microfleece fabric wicks moisture during use
- » Can be applied using only one hand
- » Allows full use of the fingers while limiting motion in the wrist

### Sizing Information

#### Circumference of Extended Wrist

Size	mm	Left Code	Right Code
XSmall/Small	100 - 150	54231L	54231R
Medium/Large	150 - 200	54233L	54233R
XLarge/XXLarge	200 - 250	54235L	54235R

**Beagle**  
ORTHOPAEDIC