


BioSkin Q Brace With Conforma Hinge



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Q Brace allows for multidirectional traction of the patella, making it a great option for anyone with patellar tracking problems or patellofemoral pain. The T-strap can be positioned on either the medial or lateral side to prevent subluxation and provide dynamic control of the patella. The conforma hinges are lightweight, flexible and provides extra support to the entire joint. Breathable and hypoallergenic materials make this patella tracking brace comfortable for all-day use. The Q Brace with conforma hinge provides excellent compression without causing discomfort. The Q Brace is available as a pull on or wrap around.

Pull on

Material - Ultima

Our pull-on knee braces are made with a highly compressive lycra-based fabric to maximise pain relief, control swelling and enhance proprioception.



Indications for use

- » Patella Tracking Disorders
- » Lateral Tracking Patella
- » Glide and Tilt Control
- » Medial Tracking Patella

Features and benefits

- » Universal design (Q strap works on either knee and can track in multiple directions)
- » Prevent subluxation
- » May also be worn as a hinged knee brace without the Q strap
- » Also available without a hinge
- » Adjustable patella tracking strap
- » Comfortable features

Sizing Information

Circumference with knee fully extended			
Size	Mid-patella	180mm above	Code
XSmall	310 - 330	410 - 460	413600
Small	330 - 360	460 - 510	413610
Medium	360 - 410	510 - 560	413620
Large	410 - 430	560 - 640	413630
XLarge	430 - 460	640 - 710	413640

Wrap Around

Material - Ultima

Many of our braces are available with our unique wraparound design. The entire front panel is Velcro compatible, making it easy to achieve a perfect fit on virtually any leg. This is also a great option for patients with limited flexibility or hand strength.



Indications for use

- » Patella Tracking Disorders
- » Lateral Tracking Patella
- » Glide and Tilt Control
- » Medial Tracking Patella

Features and benefits

- » Universal design (Q strap works on either knee and can track in multiple directions)
- » Prevent subluxation
- » May also be worn as a hinged knee brace without the Q strap
- » Also available without a hinge
- » Adjustable patella tracking strap
- » Comfortable features

Sizing Information

Circumference with knee fully extended			
Size	Mid-patella	180mm above	Code
XSmall	310 - 330	410 - 460	51380
Small	330 - 360	460 - 510	51381
Medium	360 - 410	510 - 560	51382
Large	410 - 430	560 - 640	51383
XLarge	430 - 460	640 - 710	51384