

BioSkin Q Lok Brace



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Q Lok brace helps to control patellar tracking disorders and patellofemoral pain. This knee brace comes as a pull on or wrap around and is designed to provide dynamic traction control, joint stability, and pain relief. A unique strapping system prevents patellar dislocations and incorrect tracking by applying a dynamic force to the patella, seating it correctly in the trochlear groove. This allows the joint to move throughout its range of motion smoothly and without pain. The Dynamic Traction Strap can be removed and flipped around to fit and protect either the left or right knee and prevent medial or lateral tracking disorders. Lightweight hinges on both sides of the brace anchor down the Dynamic Traction Strap and provide extra support and stability to the knee joint. These hinges naturally contour to the shape of your leg for optimized comfort. Breathable and hypoallergenic materials make the brace comfortable for all day use and provide medical-grade compression to reduce swelling, alleviate pain, and enhance recovery.

Pull on

Material - Stratus

Our pull-on knee braces are made with a highly compressive lycra-based fabric to maximize pain relief, control swelling and enhance proprioception.



Indications for use

- » Patella Tracking Disorders
- » Patella Dislocation
- » Patellofemoral Pain
- » Glide and Tilt Control

Features and benefits

- » Traction strap controls patellar glide and tilt
- » Universal left or right leg
- » Non-elastic straps control medial/lateral and anterior/posterior movement of joint
- » Rigid, aluminium uprights are flared to contour to the leg's natural shape and can be bent for additional customisation
- » SlimLok hinge: lightweight, adjustable ROM hinge with three flexion and extension stops at 0, 15 and 30

Sizing Information

Circumference with knee fully extended

Size	Mid-patella	180mm above	Code
Small	330 - 360	460 - 510	61641
Medium	360 - 410	510 - 560	61642
Large	410 - 430	560 - 640	61643
XLarge	430 - 460	640 - 710	61644

Contact customer service for other sizes.

Wrap Around

Material - Ultima

The BioSkin QLoc brace is available with our unique wraparound design. The entire front panel is hook and loop compatible, making it easy to achieve a perfect fit on virtually any leg. This is also a great option for patients with limited flexibility or hand strength.



Indications for use

- » Patella Tracking Disorders
- » Patella Dislocation
- » Patellofemoral Pain
- » Glide and Tilt Control

Features and benefits

- » Traction strap controls patellar glide and tilt
- » Fits left or right leg
- » Non-elastic straps control medial/lateral and anterior/posterior movement of joint
- » Rigid, aluminium uprights are flared to contour to the leg's natural shape and can be bent for additional customisation
- » SlimLok hinge: lightweight, adjustable ROM hinge with three flexion and extension stops at 0, 15 and 30

Sizing Information

Circumference with knee fully extended

Size	Mid - patella	180mm above	Code
Small	330 - 360	460 - 510	51641
Medium	360 - 410	510 - 560	51642
Large	410 - 430	560 - 640	51643
XLarge	430 - 460	640 - 710	51644

Contact customer service for other sizes.