

BioSkin Tennis Elbow Band


Material - Latex Free



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Tennis Elbow Band is constructed from a breathable and compressive material for comfort and effectiveness. Designed to sit comfortably just below the crease of the elbow and provide pain relieving compression to the sore area. The design includes a silicone pad and extra elastic strap for additional, targeted compression where you need it most. Position the silicone pad on the outside of the elbow to relieve the pain from lateral epicondylitis (Tennis Elbow), or move the silicone pad to the inside of the elbow for relief from medial epicondylitis (Golfer's Elbow).



Indications for use

- » Elbow pain
- » Golfer's Elbow
- » Tennis Elbow

Features and benefits

- » Constructed with BioSkin's soft, moisture-wicking materials
- » Controls joint swelling and facilitates recovery from minor injuries
- » Silicone gel pad with strap provides targeted compression to the injured area
- » Fits entirely below the crease of the elbow

Sizing Information

Circumference of largest part of forearm

Size	(mm)	Code
XSmall	200 - 230	42500
Small	230 - 260	42501
Medium	260 - 280	42502
Large	280 - 310	42503
XLarge	310 - 330	42504

Contact customer service for other sizes.

Beagle
ORTHOPAEDIC