

# BioSkin Thumb Spica


**Material** - Soft Microfleece



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

The Thumb Spica helps ease the pain and speed recovery of many wrist and thumb injuries. The thumb brace is designed to be lightweight and comfortable while providing the necessary support to the injured thumb joints. Soft microfleece is placed against the skin and a customisable, aluminum stay is used along the thumb joints for immobilization. The design allows for full use of your fingers while immobilising the injured thumb and providing compression to the wrist joint for extra blood flow and enhanced recovery. Since there are no laces, the wrist Thumb Spica is quick and easy to apply with one hand.



## Indications for use

- » Thumb pain
- » Gamekeeper's thumb
- » Skier's thumb
- » Thumb sprains

## Features and benefits

- » Customisable aluminium stay bends to fit more comfortably or reposition the thumb
- » Soft microfleece fabric wicks moisture during use
- » Can be applied using only one hand
- » Adjustable compression wrap around the wrist

### Sizing Information

#### Length of thumb

Size	mm	Code
XSmall - Medium	80 - 110	54502
Large - XXLarge	110 - 140	54505