

BioSkin TriLok Ankle Brace


Material - Ultima Latex Free



Cold Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The TriLok Ankle Brace is a slim fit design consisting of 3 components; TriLok Undersleeve, Stirrup strap and FootLok™ strap, which mimics the natural tendons in the foot to secure the ankle and allow for a safe range of motion. The TriLok is an easy, effective way to relieve the debilitating pain caused by plantar fasciitis and PTTD and also provides unrivalled protection against ankle sprains. The straps are removable for the perfect fit. Universal for left or right.



Indications for use

- » Plantar Fasciitis
- » Ankle Sprains
- » PTTD/Post-Tib Tenonitis
- » Peroneal Tendon Pain
- » Arthritis

Features and benefits

- » FootLok™ strap supports arch from above
- » Inner sleeve produces compression for swelling control and pain relief
- » Figure-8 stirrup strap provides hindfoot stabilization
- » No laces for easy application
- » Soft, antimicrobial fleece against the skin for maximum comfort
- » Designed to work as an extra ligament to improve foot and ankle problem
- » Suitable for use with socks, most shoes, without shoes and while playing sports
- » Easy to use
- » Universal left or right
- » Removable straps

Sizing Information

Size up or down for wide or narrow feet

Size	Shoe Size	Code
XSmall	2 - 3.5	53600
Small	4 - 5.5	53601
Medium	6 - 8.5	53602
Large	9 - 11.5	53603
XLarge	12 - 16	53604
XXLarge		53605
Custom		53609

Beagle
ORTHOPAEDIC