

BioSkin Visco Knee Skin **With Conforma Hinge**


Material - Latex Free



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Visco Knee Skin is an all-purpose knee brace, perfect for relieving pain from a number of injuries, including swollen joints and mild osteoarthritis pain. A silicone gel pad encircles the kneecap, providing passive massage to remove swelling and alleviate joint pain. Flexible side stays provide comfortable support to the joint. The brace is constructed from BioSkin's proprietary hypoallergenic, breathable and compressive materials meaning it's comfortable to wear all the time, while still providing medical-grade compression. Finger loops are included for easy application.



Indications for use

- » Knee instability
- » Mild knee arthritis pain
- » Patellofemoral pain

Features and benefits

- » Soft, silicone gel ring around the patella reduces swelling and increases compression around the joint
- » Also available without a hinge; contact customer service

Sizing Information

Circumference with knee fully extended

Size	Mid-patella (mm)	180mm above mid-patella (mm)	Code	With straps codes
XSmall	310 - 330	410 - 460	41700	41710
Small	330 - 360	460 - 510	41701	41711
Medium	360 - 410	510 - 560	41702	41712
Large	410 - 430	560 - 640	41703	41713
XLarge	430 - 460	640 - 710	41704	41714

Contact customer service for other sizes.

Beagle
ORTHOPAEDIC