

# Orthomerica Newport Ultra



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

The Orthomerica Newport Ultra Hip Orthosis is designed to provide protection and support for patients with problems related to hip instability. The Ultra Hip brace protects the hip joint capsule by providing a limited and safe range of motion following surgical procedures. A laterally placed abduction paddle helps to maintain the head of the femur in the acetabulum. The compression belt designed with double pull tabs anchors the hip joint and thigh cuff to promote optimal alignment. The femoral widget on the hinge offers extra adjustability for those patients who are more difficult to fit. The brace provides the user with greater mobility and support when toileting than any other hip brace. The universal fit means it can be used for the left or right and accommodates a waist size ranging from 25" – 52".



## Indications for use

- » Minimally invasive hip procedures
- » Hip resurfacing procedures
- » Labral repair, with or without gluteus medius repair
- » Patients at risk to dislocate
- » Total hip replacement
- » Post-operative or injury repair hip patients who will wear the orthosis for protection during healing
- » After repositioning of hip luxation
- » Muscular imbalance
- » Patients who have been, or are likely to be non-compliant with a rigid, bulky hip orthosis
- » Patients whose cognitive function impairs their ability to follow hip precautions
- » Arthroscopic hip procedures
- » Post-operative proximal

### Sizing Information

**Universal size fits waist 25" - 52"**

Size	Code
Universal	B/3663

## Features and benefits

- » No tools required to adjust the brace
- » Locking screw provided to prevent the patient from adjusting flexion/extension range
- » Proprioceptive reminder
- » Easily customized to provide an exact fit
- » Well-tolerated by patients
- » The virtual joint is easily set to control hip flexion, extension, adduction and abduction
- » Kinaesthetic reminder
- » Easily customized to provide an exact fit
- » Lightweight
- » Compression belt design with pull tabs anchors the hip joint and thigh cuff to promote optimal alignment
- » Universal for the right or left side
- » One size fits most
- » Can be hand washed and air dried

## Contraindications

- » Global instability
- » Patients at risk for an anterior dislocation
- » Chronic dislocation